


MOTIONMD MOBILE COMPANION APP

The Mobile Companion App was developed to add enhanced functionality and features to MotionMD including Multipage Document Upload to Patient Agreements, Inventory Counts*, Inventory Check-ins, and Inventory Transfers. The Companion App can only be used on iOS mobile devices (iPad & iPhone) at this time. There is no need to download the Companion App from the App Store from the mobile device. A prompt will show to download the App the first time the Companion App button is selected. The Companion App can be used as an adjunct when using MotionMD on the desktop via a QR code or directly opening the app from an iPad.


*Refer to the Inventory Counts knowledge guides for Companion App instructions.

Multipage Document Upload


On Patient Agreement (PA) page, under Attachments click > Add button.



In modal, click > Use Companion App button.



Select Attachment Type from the dropdown.



MOTIONMD MOBILE COMPANION APP

Multipage Document Upload (Continued)

After selecting Attachment Type
Scan the QR code to open the
Companion App on mobile.

From desktop - scan the
QR code with mobile
device from the
Companion App or from
the device camera

From mobile - after
selecting Attachment
Type, Companion App
opens automatically if it
has been downloaded

Use Companion App

Attachment Type (Required)

Chart Notes

File Attachment (Required)



Open the MotionMD
Companion App on
your iOS device. Scan
the QR code to the
left.



Download on the
App Store

Cancel Companion App

Cancel Attach File


In Companion App,
Click > Add Pages button to
Begin scanning.

After all pages have been
Scanned, click > Add X Pages
button on bottom right which
opens a Preview page.



MOTIONMD MOBILE COMPANION APP

Multipage Document Upload (Continued)



After Success! Screen, the Patient Agreement will automatically refresh and attachments are visible.

At this point, another QR code can be scanned to add additional documents, or you can navigate back to the PA page if you are on a mobile device by clicking the link at top left.



MOTIONMD MOBILE COMPANION APP

Check In Restocks

Navigate to a Restock Check In page and click > Use Companion App button.

Check In 74390 Vendor OfficeCare Standard
Restock for FMC - Brandon

Search/Scan Item Use Mobile Companion App

FMC - Brandon 5 Items Set All to Max

Product	Item #	UoM	On Hand - Outstanding (Stocking Units)	Previously Checked In (UoM)	Previously Cancelled (UoM)	Quantity Open (UoM)	Quantity to Check In (UoM)
AirSelect Short DIO - Airseal - Size S	01ES-S	10 (1 Stocking Units)	117 of 98	0	0	33	0 Max: 33
Drytex Economy Hinged Knee Sleeve DIO - Drytex - Sizes S, Small	11-0670-1	EA (1 Stocking Units)	3 of 3	0	0	2	0 Max: 2
Lace-Up Wrist Support DIO - ProCare - Size Large 7" Length - Orientation: Right	79-87297	EA	1000 of 1000	0	0	1	0 Max: 1
OA Reaction Web DIO - Drytex - Sizes Medium - Orientation: Right Medial / Left Lateral - Color: Gray	11-7426-3	EA	2 of 2	0	0	1	0 Max: 1
SRB IQ DIO - Drytex - Size Large	11-4016-4-13000	EA (1 Stocking Units)	10 of 10	0	0	9	0 Max: 9

Cancel Check In

Once you select Use Mobile Companion App – it will open up a screen with a QR code. Scan the QR code with your iPad or iPhone if you are on a computer. If you are on an iPad or iPhone, the Companion App opens automatically.

In the App's Check In screen, select either Scan Item or Manual Input.

Barcode SKUs will show as a list with the quantities scanned. Select **Upload** option when complete.

It will show a success screen and the check in counts will update on the Restock page.

10:41
Camera

Cancel Check In Upload

Scan Item
Manual Input

01EF-S 6 >

10:42
Camera

Success!

Inventory check-in was successfully submitted. You may now return to MotionMD.

Scan Another QR Code

MOTIONMD MOBILE COMPANION APP

Check In Restocks (Continued)

If a barcode was scanned or manually entered that is not a part of the Restock an error message will appear at the top of the Check In page.

Check In 74393 Vendor OfficeCare Standard

Restock for Sports Orthopedic & Spine - Escondido

Item # 11-1569-3 not found on Restock

Sports Orthopedic & Spine - Escondido 4 Items

Product	Item #	UoM	On Hand + Outstanding (Stocking Units)	Previously Checked In (UoM)	Previously Cancelled (UoM)	Quantity Open (UoM)	Quantity to Check In (UoM)
AirSelect DJO - Aircast - Size: Large	01EF-L	(1 Stocking Units)	10 of 10	0	0	1	0 Max: 1
AirSelect DJO - Aircast - Size: Medium	01EF-M	(1 Stocking Units)	50 of 50	0	0	6	0 Max: 6
AirSelect DJO - Aircast - Size: Small - Color: Chartruese	01EF-S	EA (1 Stocking Units)	5 of 5	0	0	12	6 Max: 12
Squared Toe Post-Op Shoe DJO - ProCare - Size: Small	79-81233	EA	10 of 10	0	0	4	0 Max: 4

Inventory Transfers

Inventory Transfers can be completed on the Inventory Location page by clicking the New Inventory Transfer button under the actions section on the left or by clicking Add to Transfer button on a SKU.

Sports Orthopedic & Spine - Vista

Restock 1 Items

Manual Restock

View All Restocks

View/Add Usage

Inventory Counts

New Inventory Transfer

Download Stock

Stock by Item Incomplete Restocks

01ef ☐ Needs Restock

Displaying all 4 stock records

	UoM	PAR	On Hand	Outstanding	Qty Needed
AirSelect DJO - Aircast/01EF-L Large	EA	4	3 of 4	1	0
OfficeCare Standard					
AirSelect DJO - Aircast/01EF-M	UoM	PAR	On Hand	Outstanding	Qty Needed

Add to Transfer


MOTIONMD MOBILE COMPANION APP

Inventory Transfers (Continued)

Once you open the transfer side bar, select **Use Companion App**

A modal with a QR code will open. Scan the QR code with your iPad or iPhone.




If you are already on an iPad or iPhone, the app will open automatically.



On the Transfer Screen, click **Scan Item** or **Manual Input**. Barcode SKUs will show as a list with the quantities scanned.

Select **Upload** when complete. It will show a success screen when complete and the SKUs and quantities will show in the transfer panel.


If a SKU without a stock record is uploaded an Error Message will show after upload.



MOTIONMD MOBILE COMPANION APP

Additional Notes

The Companion App can be downloaded From the App Store by searching for '**MotionMD Companion App.**'



The Companion App can also be opened by scanning the QR code with the mobile camera and clicking the link to open the app.

